



schweizerisches implantat-register  
registre suisse des implants



# Simple Export User Guide

---


V1.0 August 2020

## Simple Export by Forms

### Introduction

Simple export by forms is a feature that allows the users to run searches and export the results for the forms available.

The simple export functionality is assessable from the drop down menu as shown in the screen shot below or by clicking on the “Data Export” icon on the dashboard.



The results produced from searches are based on the location the user has logged into.

If a user logs in at clinic level, the results displayed will be those available at clinic level and the departments within the clinic. If the user is logged in at department level, they will only see results for the department they are logged into.

### Available forms

The forms available for running the searches are as follows:

#### For SIRIS Spine Users

1. Surgery SIRIS

#### For SIRIS Spine + Spine Tango users

1. Surgery SIRIS + Tango
2. Surgery 2017
3. Surgery 2017 with Trauma add-on
4. Surgery 2011
5. Surgery 2006
6. Surgery 2005
7. Staged 2011
8. Staged 2006
9. Staged 2005
10. Conservative 2018
11. Conservative 2011
12. Follow-up 2017
13. Follow-up 2011
14. Follow-up 2006

15. Follow-up 2005
16. COMI low back surgical
17. COMI neck surgical
18. COMI low back conservative
19. COMI neck conservative
20. COMI low back - baseline
21. COMI low back - 1 year post op surgery
22. COMI low back – post op
23. COMI neck - baseline
24. COMI neck - 1 year post op
25. COMI Neck – post op
26. EuroQol-3D
27. ODI
28. NDI
29. Begleiterkrankungen Patient

### Criteria for searching

To run a search for any of the forms, the key dates have to be entered in the section specified for key dates

Simple export by forms

Enter Key Date Between Filters

 and

The key dates are the date range used for performing the searches. The type of date used for searching each of the forms varies. When selecting a form to run, the key dates that apply to the form selected are as stated in the table below.

### SIRIS Spine

	Form Name	Key Date
1	Surgery SIRIS	Surgery Date

### SIRIS Spine and Tango

Nr	Form Name	Key Date
1	Surgery SIRIS + Tango	Surgery Date
2	Surgery 2017	Surgery Date
3	Surgery 2017 with Trauma add-on	Surgery Date
4	Surgery 2011	Surgery Date
5	Surgery 2006	Surgery Date
6	Surgery 2005	Surgery Date
7	Staged 2011	Surgery Date

8	Staged 2006	Surgery Date
9	Staged 2005	Surgery Date
10	Conservative 2018	Therapy Date
11	Conservative 2011	Therapy Date
12	Follow-up 2017	Follow Up Date
13	Follow-up 2011	Follow Up Date
14	Follow-up 2006	Follow Up Date
15	Follow-up 2005	Follow Up Date
16	COMI low back surgical	Evaluation Date
17	COMI neck surgical	Evaluation Date
18	COMI low back conservative	Evaluation Date
19	COMI neck conservative	Evaluation Date
20	COMI low back - baseline	Evaluation Date
21	COMI low back - 1 year post op surgery	Evaluation Date
22	COMI low back – post op	Evaluation Date
23	COMI neck - baseline	Evaluation Date
24	COMI neck - 1 year post op	Evaluation Date
25	COMI neck – post op	Evaluation Date
26	EuroQol-3D	Examination Date
27	ODI	Examination Date
28	NDI	Examination Date
29	Begleiterkrankungen Patient	Untersuchungsdatum

Once a date range has been selected, click on the ‘run button’ (circled below) for the required form to execute the search.

Summary	Searches
 	1. Surgery 2017 - Filter by Surgery date
 	2. COMI low back surgical - Filter by Date of evaluation
 	3. COMI neck surgical - Filter by Date of evaluation
 	4. COMI low back conservative - Filter by Date of evaluation

Once the search has been completed, a number will be displayed next to the form with the results found.

Summary	Searches
 6 	1. Surgery 2017 - Filter by Surgery date

Enter Key Date Between Filters

01/01/2019  and 01/01/2020 

To re-run your search for another form using the same date range, click the run button next to the required form and the results will be displayed.

**Simple export by forms**

Enter Key Date Between Filters

  and  

Summary			Searches
	6		1. Surgery 2017 - Filter by Surgery date
	5		2. COMI low back surgical - Filter by Date of evaluation

Where there are no results, the number 'zero' will be displayed.

**Simple export by forms**

Enter Key Date Between Filters

  and  

Summary			Searches
	6		1. Surgery 2017 - Filter by Surgery date
	5		2. COMI low back surgical - Filter by Date of evaluation
	0		3. COMI neck surgical - Filter by Date of evaluation

If the date range is changed, the results will be cleared and the searches will have to be re-run with the new dates. The results shown are always based on the key dates specified when the search is run.

## Exporting Results

To export the results of any form that has been searched, click on the export button next to the number of results.

**Simple export by forms**

Enter Key Date Between Filters

  and  

Summary			Searches
	6		1. Surgery 2017 - Filter by Surgery date
	5		2. COMI low back surgical - Filter by Date of evaluation
	0		3. COMI neck surgical - Filter by Date of evaluation

A zipped file will be created with the results of your search and downloaded to your device.

